**memosplint**<sup>®</sup>  
The first thermoelastic blank for Dental

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deutsch | english | français | italiano | español

Stand: 01/2021

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Content | Index | Indice | Indice

deutsch .....	Seite 4 – 7
english .....	Page 8 – 11
français .....	Page 12 – 15
italiano .....	Página 16 – 19
español .....	Pagina 20 – 23

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Phone: +49 (0) 231 - 55 56 - 0

Fax: +49 (0) 231 - 55 56 - 30

## **Thermoplastic bite guard splints from optimill memosplint**

The bite guard splints milled from optimill memosplint show thermoplastic flexibility. The result for patients is exceptional stress-free wearing comfort thanks to maximum adaptation to the tooth situation. What's more, the transparent flexible bite guard splints are highly resistant to breakage and highly transparent for excellent visibility.

### **Indications for optimill memosplint**

Dental resins for the production of flexible dental therapeutic occlusal splints of all kinds

### **Contraindications for optimill memosplint**

The material should not be used for any purpose other than the generative CAD/CAM fabrication of dental bite guard splints. Optimill memosplint is not suitable for fabricating denture bases. Undesirable biological reactions (e.g. allergies to material components) may occur in very rare cases. Do not use bite guard splint made with optimill memosplint if incompatibilities or allergies to the material components are known.


### **Wall thicknesses / CAD temporary cast**

When splints are being designed, the wall thickness must not be less than the minimum of 1.0 mm. Please check this requirement in your design software. In the occlusal region, the thickness of the splints after trimming must not be less than 0.9 mm. A reduction to 0.9 mm in the labial region is permissible for aesthetic reasons. Any corners and edges must be rounded off.

### **Milling preparation**

During nesting, please ensure that enough attachments are used. (It is a thermoplastic material and would otherwise become deformed in the machine. This would result in the burs being broken).

## Milling

A conventional dry milling strategy for PMMA can be used. Single-edged PMMA tools should be used in the machine. The attachments are cleaned as usual after milling.

## Post-processing of objects

The splints fabricated from optimill memosplint can be processed with suitable cross-toothed carbide cutters or suitable cutting wheels for resin. Final polishing of the splints is essential to prevent any accumulation of plaque on them. Polishing is done with standard polishing brushes and polishing pastes. Please be aware that excessive contact pressure during polishing may overheat the material and cause burns.

## Cleaning the objects

Before further processing, the objects should not be cleaned with excessive evaporation due to the high heat generation. An ultrasonic bath with a maximum temperature of about 35°C and a soft toothbrush are more suitable. Cleaning solutions and denture cleaners are not recommended.

## Technical Data

Chemical composition in %	
Poly(m)ethylacrylate	> 90
Cross-linked copolymers of methacrylic acid	> 90
1,2-cyclohexane dicarboxylic acid diisononyl ester	> 10
Material data	
Flexural strength (23°C)	> 60 MPa
Flexural strength (37°C)	< 60 MPa
Density	approx. 1.1-1.2 g/cm <sup>3</sup>
Color	clear transparent

## Notes

Restorations are custom-made products according to Directive 93/42/EEC. Store in clear water at least 12 hours, preferably overnight before delivery. For each process, indicate the LOT number required to identify the material.

## Storage

Always store in clear water, when not in use. Avoid direct sunlight. No further special measures required.

## Hazard warnings

Processing optimill memosplint produces dust that may irritate the eyes, skin and respiratory tract. So please always ensure that the extractor fan at your workplace is working flawlessly. See safety data sheet on [www.dentona.de](http://www.dentona.de).

## Disposal

Restorations made with optimill memosplint are insoluble in water, inactive, pose no threat to ground water and can thus be disposed of in normal domestic waste.

## Warranty

Our application-specific recommendations, no matter whether they are given orally, in writing or by way of practical tutorials, are based on our own experiences and tests and may therefore only be considered as recommended values. Our products are subject to continuous enhancements. We therefore reserve the right to make changes in design and composition.

## Use

Only to be used by professionals. Rx only.

Patient population: adults and adolescents

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dentona AG
**dentona AG** | Otto-Hahn-Str. 27 | 44227 Dortmund